

Intraoperative Ophthalmic Checklist	√ when completed
TIME OUT!	
RECORD Color or number of each instrument and tray used	
Use <u>disposable</u> intraocular cannula's	
Intraocular medications should be preservative-free.	
Utilize the yellow Viscoat "cannula locking device"	
Label all syringes and basins	
Do not draw BSS w/ EPI from the bag on the Infiniti tray	
Verify IOL power w/ surgeon	
Do not touch IOL with the gloved hand	
Do not touch any instrument tip with a gloved hand	
Wipe instruments and flush or immerse in sterile distilled water in the O.R. immediately following use	
Flush Ozil handpiece with 120ml sterile distilled water utilizing the lumen adapter.	
Flush Ozil lumens with 60ml of air prior to sterilization.	
Extend plunger on Monarch IOL inserter fully and wipe w/ Weck wipe	
Flush Monarch IOL inserter with sterile distilled water and leave the instrument to soak in water	
Keep all used ophthalmic instruments moist with sterile water	
Recommendation: Refrain from using ophthalmic ointments and all ophthalmic medications with preservatives	